



MediRay Healthcare (P) Limited

Committed to Customer Committed to the future...**MediRay**

Ph: +91-11-47063883, 27186595, Fax: +91-11-47091312

Web: www.mediray.in, Mail: info@mediray.in,

Skype: MediRayIndia, MSN: MediRay@live.com

UltimaX- 50

Specification Of X-Ray Tube Head

Output	50mA at 50 kV 40mA at 65kV 35mA at 65 kV 30mA at 85 kV
Time	Microprocessor 0.1 to 6 sec. 60 steps
Tube	1.8mm focus Bel dsa-1or Equivalent
Input Compensation	230V. 5AMP. 210To 250V.
Beamlimita-	Light Beam Collimator
Weight	20 Kgs.



Special Features

- ☒ Truly Portable Carry away box type Machine
- ☒ Voltage Over Load Protection
- ☒ Tape built FFD Measurement
- ☒ Micro Processor Based Digital with Accuracy Of 0.01Sec.
- ☒ Exposure Time Can be set by soft buttone
- ☒ Full P2oof Radiation Protection.
- ☒ Line Voltege Compensation
- ☒ Hand Operatad switch is provided for X- Ray exposure in single touch

STANDARD ACCESSORIES:

Power Supply Cable, Exposure Switch with Cable Carrying Case Beam limitation Device

STANDS

ALLO Ratche Stand:- Light Wight (12Kg.) Dismantable Stand Stardy Convenient Ratch Typestand.